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Cat:C100603

Size:2x50mL/2x100mL

Storage:2- 8°C, avoid light, valid for 6 months.

Kit Components:

Reagent		2x 50mL	2x 100mL	Storage
Reagent(A):ModifiedOil Red O Stain Solution	A1: Oil Red O Stain A	30mL	60mL	2-8°C,avoid light
	A2: Oil Red O Stain B	20mL	40mL	RT
Mix A1 and A2 with the ratio of 3:2 and place for 10 min to form Modified Oil Red O Stain Solution.Configure before use.				
Reagent(B):Mayer's Hematoxylin Solution		50mL	100mL	2-8°C,avoid light

Introduction:

Neutral fat stains often use Sudan II, Sudan III, Sudan IV, Sudan black B, oil red O and so on. Sudan stain is often used in traditional methods. Recently, azo dye oil red O is more suitable for dyeing fat. Oil red O is a strong lipid solvent and dye. It is easy to bind to triglycerides in droplet shape, but has a slightly poor binding to phospholipids.

Modified Oil Red O Stain Kit is mainly used to show fatty degeneration of tissues and organs and abnormal lipid-like sedimentation. It often occurs in fatty degeneration of liver, kidney, heart and other parenchymal organs that most of the neutral fat droplets appear in cells. It can identify and diagnose tumors and their properties in adipose tissue. The positive staining results of fat were orange to red, but the specific color depended on the concentration of lipid.

Self Provided Materials:

60% isopropanol, distilled water, bluing solution, glycerol gelatin

Protocol:(for reference only)

1. Cut frozen sections 10 um thick,unfix or wash by water after fix in 10% formalin for 10min.Rinse sections slightly with distilled water.
2. Soak sections in 60% isopropanol for 20-30s.
3. Stain section in Modified Oil Red O Stain Solution (capped) for 10-15 min.
4. Colour separation: Wash slightly for 5-10s in 60% isopropanol to remove the dye solution. Rinse sections with distilled water for 1 min.

5. Re-dyeing by Mayer's Hematoxylin Solution for 1-2 min.
6. Rinse by tap water for 10 min or blue with bluing solution for 3min.
7. Rinse sections slightly with distilled water. Absorb the surrounding water with filter paper.
8. Glycerol gelatin or arabinose glue scaling is solid.

Result:

Neutral Fat	Orange or Red
Nucleus	Blue

Note:

1. The Modified Oil Red O Stain Solution is not stable enough and easy to precipitate, so it is not suitable to prepare it in advance. Filter before use.
2. If 60% isopropanol is not available, 70% ethanol can also be used.
3. Lipids are easily soluble in organic solvents, and it is not recommended to prepare paraffin or plastic coated sections for staining. It is usually recommended to prepare frozen sections or carbon wax sections.
4. Adipose cells are usually larger than normal cells, and frozen sections should not be too thin (10-20um) Thin slices often cause lipid loss.
5. It is normal for reagent A2 to become turbid during low-temperature storage and transportation, and it can be restored to clarity and normal use by boiling water for 5 min.
6. The re-staining time of Mayer's hematoxylin staining solution should not be too long, and the staining results should not be stored for a long time. Observation and photography should be done as soon as possible.
7. If long-term preservation is required, the edge between the cover glass and the slide can be sealed with neutral resin.
8. For your safety and health, please wear laboratory clothes and disposable gloves when operating.